



**MOBILITY AIDS**

**POWER WHEELCHAIR**



## **THUNDER POWER WHEELCHAIR PW-01**

- MOTOR 180\*2 BRUSH
- BATTERY 24V 12AH LEAD ACID
- MAX LOADING 120KG
- CHARGING TIME 5-7H
- DRIVING DISTANCE 15~20KM
- N.W.(WITH BATTERY) 31KG
- N.W.(WITHOUT BATTERY) 25KG
- PACKING SIZE 68\*35\*73CM